

Trouble hearing?

Having trouble hearing in some situations?



When talking one-on-one

Do you find it hard to follow a conversation when chatting with a friend or relative?



When dining or in meetings

Do you have trouble hearing what is being said at family dinners or get-togethers?



When watching television

Do you miss words when you watch TV at normal volume?



When driving

Do you find it hard to hear in a noisy environment like for instance in your car?

This is the solution



Superior sound

Do you remember the old cassette tapes? Along with the music came unwanted noise and complicated handling. The sound technology in Maxi, however, is fully digital and uses the same technology as found in modern CD/DVD players. As a result, Maxi has crystal clear sound and can be used for hours without you growing tired of it.

Easy to use

Maxi is exceptionally user friendly. Its shape fits nicely on the table or in your hand, and the rubberized material makes it easy to grip. The controls are simple to use, even for those with impaired vision or dexterity problems.

High quality

Maxi is designed by Bellman & Symfon of Sweden and comes with a two-year warranty.

Discover Bellman® Audio Maxi



The difference.

Discover how Maxi can improve your hearing in everyday life.



How Maxi works

Maxi's built-in microphone picks up and amplifies sounds using digital technology. You hear the sound through your headphones, earphones or a neck loop if you are using hearing aids.

Maxi also has a built-in telecoil that works with inductive loop systems that are common in public areas.



When talking one-on-one



When dining or in meetings



Cable Kit - for Hi-fi TV sound

Use the cable kit to connect your Maxi to the TV and enjoy high quality sound in full stereo.



If you push the M-button while listening, you can also hear people talking without removing your listening device.

Tie-Clip Microphone

For meetings and easy TV-listening

Specifically developed for Maxi as an easy-to-setup hearing solution in mono. Just place the microphone close to the TV's loudspeaker. The Tie-Clip Microphone can also be used at dinners or meetings to improve hearing even further.



When driving



Listening alternatives



Headphones - for easy use

Stereo headphones with superior sound quality. The headphones are easy to use, foldable and very comfortable to wear.



Earphones - for a discrete look

Stereo earphones with excellent sound. The earphones are light-weight, discrete and delivers sound with high clarity and great detail.




Neck loop - for your hearing aid



The neck loop sends sound to hearing aids equipped with a telecoil. You can also place Maxi close to a speaker or connect it to the TV with the cable kit. Maxi will function as a remote amplified microphone and make the sound much clearer.

When watching television



A car is one of the hardest environments to conduct a conversation in. With Maxi the problem is history. Simply attach the tie-clip microphone on the person you want to talk to and enjoy every great conversation.

Maxi comes with

- Maxi Digital Communication Aid
- Listening device of your choice (headphones, earphones or neck loop)
- Audio cable kit for TV
- Two-year warranty