

# cm-BT2




## **The Bluetooth headset CM-BT2:**

**Keep your hands free when using your mobile phone, receive music from your PC directly at your ear, always in your preferred volume ...**

**... the CM-BT2 is your wireless audio connection to many Bluetooth-enabled devices.**

HUMANTECHNIK –  
Your partner for audiological products





Want to listen to music in stereo from your PC? Have your navigation system speak the route directly into your hearing aid? Enjoy even greater convenience and understand everything better on your mobile phone? **The CM-BT2: Modern Bluetooth application plus intelligent audio-logical technology for relaxed hearing and excellent understanding.**



The Bluetooth headset, CM-BT2 (1:1)

As soon as an incoming call is detected, your Bluetooth mobile phone sets up a radio link to your CM-BT2. With the touch of a single button, you can answer all your calls via the audio amplifier and immediately have your hands free once again. You can speak easily and understand your calling partner better than you ever have before.



The CM-BT2 set includes the Bluetooth hearing amplifier with teleloop and rechargeable battery, two chargers for 230 V and 12 V and a USB-charging cable.

The following **accessories** are available: Stereo earphones, light-weight headphones, induction link (monaural / binaural) audio cable (monaural / binaural).

#### What is Bluetooth?

»Bluetooth« is an internationally standardized data interface for radio applications over short distances. Many manufacturers of consumer electronics, communications technology and computers now offer products equipped with Bluetooth interfaces. Bluetooth technology is most common today in combination with the hands-free operation of mobile telephones.

#### Complete Bluetooth comfort and more!

The Bluetooth headset CM-BT2 is specially developed for the hearing impaired. It allows infinitely variable

volume control and offers a differentiated range of sound with the treble/bass switch.

The CM-BT2 uses the teleloop to transmit audio signals inductively to hearing aids equipped with a telephone coil. In addition, the amplifier also provides a connection jack for stereo ear- and headphones. (Noise level under earphone use up to 118 dB).

#### Versatile and convenient utilization

Apart from being a mobile phone hands-free system and its audio-amplifying abilities, the CM-BT2 also has an inbuilt microphone.

Therefore it enables sound to be transmitted via other Bluetooth-enabled media, for example a PC. At the same time it also offers

all the comfort of a wireless sound system for your stereo equipment or your computer.

If your PC does not have a Bluetooth interface, our



USB-Bluetooth adaptor (special accessory) provides the necessary contact to the CM-BT2. Compatible adaptors for stereo systems and more are available at your specialised dealer.



**HUMANTECHNIK**

#### HUMANTECHNIK GmbH

Im Wörth 25 · D-79576 Weil am Rhein  
 Telefon: +49 (0) 7621 / 9 56 89 - 0  
 Fax: +49 (0) 76 21 / 9 56 89 - 70  
 Internet: www.humantechnik.com  
 E-Mail: info@humantechnik.com

Your hearing aid dispenser will advise you comprehensively and competently